

Run Info

Experiment Name	ncov_ucdh_env1_run1
Sample ID	ncov_ucdh_env1_run1
Run ID	25240491-27d6-4219-a825-722520587610
Flow Cell Id	FAN33832
Start Time	May 23, 00:32
Run Length	15h 13m

Run Summary

Reads Generated	405.97 K
Bases Generated	214.86 Mb
Estimated Bases	225.11 Mb


Run Parameters

Flow Cell Type	FLO-MIN106
Kit	SQK-LSK109
Basecalling	on
Specified Run Length	72 hours
Initial Bias Voltage	-180 mV
FAST5 Output	Enabled
FAST5 Output Options	zlib_compress,fastq,raw
FAST5 Reads per File	1000
FASTQ Output	Enabled
FASTQ Reads per File	1000
Active Channel Selection	Enabled
Mux Scan Period	1 hour 30 minutes
Reserved Pores	0 %
Basecall Model	Fast basecalling


Versions

MinKNOW Core	3.6.5
Bream	4.3.16
Guppy	3.2.10

Cumulative Output Reads




Cumulative Output Bases




Read Length Histogram Estimated Bases

Estimated N50: 545 b




Read Length Histogram Basecalled Bases


Estimated N50: 510 b




Duty Time Grouped




Duty time Categorised




Mux Scan Grouped




Mux Scan Categorised




Temperature History




Bias Voltage History




Translocation Speed



QScore



Disk Write Performance



Run Debug Messages

- Flow cell FAN33832 has 1508 pores available for sequencing. Starting sequencing with 510 pores May 23, 15:39
- Performing Mux Scan May 23, 15:37
- Flow cell FAN33832 has 1511 pores available for sequencing. Starting sequencing with 492 pores May 23, 14:23
- Performing Mux Scan May 23, 14:21
- Flow cell FAN33832 has 1550 pores available for sequencing. Starting sequencing with 512 pores May 23, 12:51
- Performing Mux Scan May 23, 12:49
- Flow cell FAN33832 has 1558 pores available for sequencing. Starting sequencing with 507 pores May 23, 11:19
- Performing Mux Scan May 23, 11:17
- Flow cell FAN33832 has 1566 pores available for sequencing. Starting sequencing with 510 pores May 23, 09:47
- Performing Mux Scan May 23, 09:45
- Flow cell FAN33832 has 1587 pores available for sequencing. Starting sequencing with 512 pores May 23, 08:15
- Performing Mux Scan May 23, 08:13
- Flow cell FAN33832 has 1597 pores available for sequencing. Starting sequencing with 511 pores May 23, 06:43
- Performing Mux Scan May 23, 06:41
- Flow cell FAN33832 has 1607 pores available for sequencing. Starting sequencing with 512 pores May 23, 05:11
- Performing Mux Scan May 23, 05:09
- Flow cell FAN33832 has 1632 pores available for sequencing. Starting sequencing with 512 pores May 23, 03:39
- Performing Mux Scan May 23, 03:37
- Flow cell FAN33832 has 1659 pores available for sequencing. Starting sequencing with 512 pores May 23, 02:07
- Performing Mux Scan May 23, 02:05
- Flow cell FAN33832 has 1673 pores available for sequencing. Starting sequencing with 512 pores May 23, 00:35
- Performing Mux Scan May 23, 00:33
- Starting sequencing procedure May 23, 00:33
- Waiting up to 300 seconds for temperature to stabilise at 34.0°C May 23, 00:32
- Disk / has 1881 GB space remaining May 23, 00:32